

# UQCRC2 Conjugated Antibody

Catalog No: #C56271

Package Size: #C56271-AF350 100ul #C56271-AF405 100ul #C56271-AF488 100ul #C56271-AF555 100ul #C56271-AF594 100ul #C56271-AF647 100ul #C56271-AF680 100ul #C56271-AF750 100ul #C56271-Biotin 100ul #C56271-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	UQCRC2 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	UQCRC2 Antibody detects endogenous levels of total UQCRC2
Immunogen Description	A synthesized peptide derived from human UQCRC2
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Core protein II; QCR2; mitochondrial; Ubiquinol cytochrome c reductase core protein II; UQCR2; Uqcr2;
Accession No.	Uniprot:P22695
Calculated MW	48kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

## Product Description

This is a component of the ubiquinol-cytochrome c reductase complex (complex III or cytochrome b-c1 complex), which is part of the mitochondrial respiratory chain. The core protein 2 is required for the assembly of the complex.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.