

TRADD Conjugated Antibody

Catalog No: #C56274



Package Size: #C56274-AF350 100ul #C56274-AF405 100ul #C56274-AF488 100ul #C56274-AF555 100ul #C56274-AF594 100ul #C56274-AF647 100ul #C56274-AF680 100ul #C56274-AF750 100ul #C56274-Biotin 100ul #C56274-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	TRADD Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	TRADD Antibody detects endogenous levels of total TRADD
Immunogen Description	A synthesized peptide derived from human TRADD
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	AA930854; TNFR1 associated DEATH domain protein; tradd;
Accession No.	Uniprot:Q15628
Calculated MW	34kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Adapter molecule for TNFRSF1A/TNFR1 that specifically associates with the cytoplasmic domain of activated TNFRSF1A/TNFR1 mediating its interaction with FADD. Overexpression of TRADD leads to two major TNF-induced responses, apoptosis and activation of NF-kappa-B.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.