

CTLA4 (CD152) Conjugated Antibody

Catalog No: #C56275

Package Size: #C56275-AF350 100ul #C56275-AF405 100ul #C56275-AF488 100ul #C56275-AF555 100ul #C56275-AF594 100ul #C56275-AF647 100ul #C56275-AF680 100ul #C56275-AF750 100ul #C56275-Biotin 100ul #C56275-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	CTLA4 (CD152) Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	CTLA4 (CD152) Antibody detects endogenous levels of total CTLA4 (CD152)
Immunogen Description	A synthesized peptide derived from human CTLA4 (CD152)
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CTLA4;CD;CD152;CELIAC3;CTLA-4;GRD4;GSE;ICOS;IDDM12;
Accession No.	Uniprot:P16410
Calculated MW	43kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Inhibitory receptor acting as a major negative regulator of T-cell responses. The affinity of CTLA4 for its natural B7 family ligands, CD80 and CD86, is considerably stronger than the affinity of their cognate stimulatory coreceptor CD28.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.