

ROC1 Conjugated Antibody

Catalog No: #C56277

Package Size: #C56277-AF350 100ul #C56277-AF405 100ul #C56277-AF488 100ul #C56277-AF555 100ul #C56277-AF594 100ul #C56277-AF647 100ul #C56277-AF680 100ul #C56277-AF750 100ul #C56277-Biotin 100ul #C56277-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	ROC1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	ROC1 Antibody detects endogenous levels of total ROC1
Immunogen Description	A synthesized peptide derived from human ROC1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	BA554C12.1; MGC13357; MGC1481; RNF75; ROC1;
Accession No.	Uniprot:P62877
Calculated MW	12kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

E3 ubiquitin ligase component of multiple cullin-RING-based E3 ubiquitin-protein ligase complexes which mediate the ubiquitination and subsequent proteasomal degradation of target proteins, including proteins involved in cell cycle progression, signal transduction, transcription and transcription-coupled nucleotide excision repair.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.