

TAX1BP3 Conjugated Antibody

Catalog No: #C56281

Package Size: #C56281-AF350 100ul #C56281-AF405 100ul #C56281-AF488 100ul #C56281-AF555 100ul #C56281-AF594 100ul #C56281-AF647 100ul #C56281-AF680 100ul #C56281-AF750 100ul #C56281-Biotin 100ul #C56281-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	TAX1BP3 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	TAX1BP3 Antibody detects endogenous levels of total TAX1BP3
Immunogen Description	A synthesized peptide derived from human TAX1BP3
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Glutaminase interacting protein 3; Tax interaction protein 1; Tax1 binding protein 3; Tax1bp3; TIP1;
Accession No.	Uniprot:O14907
Calculated MW	14kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

May play a role in the Rho signaling pathway. May act as an inhibitor of the Wnt signaling pathway. May play a role in activation of CDC42 by the viral protein HPV16 E6.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.