

# Factor X Conjugated Antibody

Catalog No: #C56283

Package Size: #C56283-AF350 100ul #C56283-AF405 100ul #C56283-AF488 100ul #C56283-AF555 100ul #C56283-AF594 100ul #C56283-AF647 100ul #C56283-AF680 100ul #C56283-AF750 100ul #C56283-Biotin 100ul #C56283-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Factor X Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	Factor X Antibody detects endogenous levels of total Factor X
Immunogen Description	A synthesized peptide derived from human Factor X
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Activated factor Xa heavy chain; Coagulation factor X; F10; factor Xa; FXA; Prothrombinase; Stuart factor;
Accession No.	Uniprot:P00742
Calculated MW	74kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

## Product Description

Factor Xa is a vitamin K-dependent glycoprotein that converts prothrombin to thrombin in the presence of factor Va, calcium and phospholipid during blood clotting.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.