

Claudin 4 Conjugated Antibody

Catalog No: #C56285

Package Size: #C56285-AF350 100ul #C56285-AF405 100ul #C56285-AF488 100ul #C56285-AF555 100ul #C56285-AF594 100ul #C56285-AF647 100ul #C56285-AF680 100ul #C56285-AF750 100ul #C56285-Biotin 100ul #C56285-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Claudin 4 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	Claudin 4 Antibody detects endogenous levels of total Claudin 4
Immunogen Description	A synthesized peptide derived from human Claudin 4
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CLD4; CLDN4; CPE- receptor; CPE-R; CPER; CPETR1; Clostridium perfringens enterotoxin receptor; WBSCR8; Williams-Beuren syndrome chromosome region 8 protein;
Accession No.	Uniprot:O14493
Calculated MW	18kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Plays a major role in tight junction-specific obliteration of the intercellular space.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.