

CD97 Conjugated Antibody

Catalog No: #C56288



Package Size: #C56288-AF350 100ul #C56288-AF405 100ul #C56288-AF488 100ul #C56288-AF555 100ul #C56288-AF594 100ul #C56288-AF647 100ul #C56288-AF680 100ul #C56288-AF750 100ul #C56288-Biotin 100ul #C56288-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	CD97 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse
Specificity	CD97 Antibody detects endogenous levels of total CD97
Immunogen Description	A synthesized peptide derived from human CD97
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	ADGRE5; Adhesion G protein-coupled receptor E5; CD97 antigen; Leukocyte antigen CD97; Seven span transmembrane protein; TM7LN1;
Accession No.	Uniprot:P48960
Calculated MW	92kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF: 1:50-1:200

FC: 1:50-1:200

Product Description

Receptor potentially involved in both adhesion and signaling processes early after leukocyte activation. Plays an essential role in leukocyte migration.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.