

BOB1 Conjugated Antibody

Catalog No: #C56290

Package Size: #C56290-AF350 100ul #C56290-AF405 100ul #C56290-AF488 100ul #C56290-AF555 100ul #C56290-AF594 100ul #C56290-AF647 100ul #C56290-AF680 100ul #C56290-AF750 100ul #C56290-Biotin 100ul #C56290-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	BOB1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	BOB1 Antibody detects endogenous levels of total BOB1
Immunogen Description	A synthesized peptide derived from human BOB1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	BOB 1; BOB-1; OBF 1; OBF1; OCA B; OCA-B; OCAB; OCT binding factor 1; Pou2af1;
Accession No.	Uniprot:Q16633
Calculated MW	35kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Transcriptional coactivator that specifically associates with either OCT1 or OCT2. It boosts the OCT1 mediated promoter activity and to a lesser extent, that of OCT2. It has no intrinsic DNA-binding activity. It recognizes the POU domains of OCT1 and OCT2. It is essential for the response of B-cells to antigens and required for the formation of germinal centers.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.