

GRB2 Conjugated Antibody

Catalog No: #C56294

Package Size: #C56294-AF350 100ul #C56294-AF405 100ul #C56294-AF488 100ul #C56294-AF555 100ul #C56294-AF594 100ul #C56294-AF647 100ul #C56294-AF680 100ul #C56294-AF750 100ul #C56294-Biotin 100ul #C56294-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	GRB2 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	GRB2 Antibody detects endogenous levels of total GRB2
Immunogen Description	A synthesized peptide derived from human GRB2
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	ASH; Grb2; Grb3 3; HT027; NCKAP2; Protein ash; SEM5;
Accession No.	Uniprot:P62993
Calculated MW	25kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Adapter protein that provides a critical link between cell surface growth factor receptors and the Ras signaling pathway.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.