

CD33 Conjugated Antibody

Catalog No: #C56297



Package Size: #C56297-AF350 100ul #C56297-AF405 100ul #C56297-AF488 100ul #C56297-AF555 100ul #C56297-AF594 100ul #C56297-AF647 100ul #C56297-AF680 100ul #C56297-AF750 100ul #C56297-Biotin 100ul #C56297-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	CD33 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	CD33 Antibody detects endogenous levels of total CD33
Immunogen Description	A synthesized peptide derived from human CD33
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CD33; gp67; My9; p67; SIGLEC3;
Accession No.	Uniprot:P20138
Calculated MW	65-75kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Putative adhesion molecule of myelomonocytic-derived cells that mediates sialic-acid dependent binding to cells. Preferentially binds to alpha-2,6-linked sialic acid.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.