

# Egr1 Conjugated Antibody

Catalog No: #C56298

Package Size: #C56298-AF350 100ul #C56298-AF405 100ul #C56298-AF488 100ul #C56298-AF555 100ul #C56298-AF594 100ul #C56298-AF647 100ul #C56298-AF680 100ul #C56298-AF750 100ul #C56298-Biotin 100ul #C56298-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Egr1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	Egr1 Antibody detects endogenous levels of total Egr1
Immunogen Description	A synthesized peptide derived from human Egr1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	TIS8; AT225; G0S30; NGFI-A; ZNF225; KROX-24; ZIF-268; EGR1;
Accession No.	Uniprot:P18146
Calculated MW	80kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

## Product Description

Transcriptional regulator. Recognizes and binds to the DNA sequence 5'-CGCCCCCGC-3'(EGR-site). Activates the transcription of target genes whose products are required for mitogenesis and differentiation.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.