

ZEB1 Conjugated Antibody

Catalog No: #C56299



Package Size: #C56299-AF350 100ul #C56299-AF405 100ul #C56299-AF488 100ul #C56299-AF555 100ul #C56299-AF594 100ul #C56299-AF647 100ul #C56299-AF680 100ul #C56299-AF750 100ul #C56299-Biotin 100ul #C56299-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	ZEB1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human
Specificity	AREB6 Antibody detects endogenous levels of total AREB6
Immunogen Description	A synthesized peptide derived from human AREB6
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	BZP; TCF8; AREB6; FECD6; NIL2A; PPCD3; ZFHEP; ZFHX1A; DELTAEF1;
Accession No.	Uniprot:P37275
Calculated MW	200kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Inhibits interleukin-2 (IL-2) gene expression. May be responsible for transcriptional repression of the IL-2 gene. Enhances or represses the promoter activity of the ATP1A1 gene depending on the quantity of cDNA and on the cell type.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.