

CD23 Conjugated Antibody

Catalog No: #C56300

Package Size: #C56300-AF350 100ul #C56300-AF405 100ul #C56300-AF488 100ul #C56300-AF555 100ul #C56300-AF594 100ul #C56300-AF647 100ul #C56300-AF680 100ul #C56300-AF750 100ul #C56300-Biotin 100ul #C56300-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	CD23 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human
Specificity	CD23 Antibody detects endogenous levels of total CD23
Immunogen Description	A synthesized peptide derived from human CD23
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CD23; FCE2; CD23A; IGEBF; CLEC4J; FCER2;
Accession No.	Uniprot:P06734
Calculated MW	36kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

This receptor has essential roles in the regulation of IgE production and in the differentiation of B-cells (it is a B-cell-specific antigen).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.