

# CD30 Conjugated Antibody

Catalog No: #C56302

Package Size: #C56302-AF350 100ul #C56302-AF405 100ul #C56302-AF488 100ul #C56302-AF555 100ul #C56302-AF594 100ul #C56302-AF647 100ul #C56302-AF680 100ul #C56302-AF750 100ul #C56302-Biotin 100ul #C56302-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CD30 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	CD30 Antibody detects endogenous levels of total CD30
Immunogen Description	A synthesized peptide derived from human CD30
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CD30; D1S166E; KI1; Tnfrsf8;
Accession No.	Uniprot:P28908
Calculated MW	120kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

## Product Description

Receptor for TNFSF8/CD30L. May play a role in the regulation of cellular growth and transformation of activated lymphoblasts. Regulates gene expression through activation of NF-kappa-B.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.