

BAG2 Conjugated Antibody

Catalog No: #C56306

Package Size: #C56306-AF350 100ul #C56306-AF405 100ul #C56306-AF488 100ul #C56306-AF555 100ul #C56306-AF594 100ul #C56306-AF647 100ul #C56306-AF680 100ul #C56306-AF750 100ul #C56306-Biotin 100ul #C56306-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	BAG2 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	BAG2 Antibody detects endogenous levels of total BAG2
Immunogen Description	A synthesized peptide derived from human BAG2
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	BAG 2; Bag2;
Accession No.	Uniprot:O95816
Calculated MW	24kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Inhibits the chaperone activity of HSP70/HSC70 by promoting substrate release.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.