

GST Conjugated Antibody

Catalog No: #C56319

Package Size: #C56319-AF350 100ul #C56319-AF405 100ul #C56319-AF488 100ul #C56319-AF555 100ul #C56319-AF594 100ul #C56319-AF647 100ul #C56319-AF680 100ul #C56319-AF750 100ul #C56319-Biotin 100ul #C56319-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	GST Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	All
Specificity	GST Antibody detects endogenous levels of C-terminal, internal, and N-terminal GST-tagged proteins.
Immunogen Description	A synthesized peptide derived from GST
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	N/A
Accession No.	Uniprot:N/A
Calculated MW	26kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

GST-Tag antibody is part of the Tag series of antibodies, the best quality in the research.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.