

GRK3 Conjugated Antibody

Catalog No: #C56322

Package Size: #C56322-AF350 100ul #C56322-AF405 100ul #C56322-AF488 100ul #C56322-AF555 100ul #C56322-AF594 100ul #C56322-AF647 100ul #C56322-AF680 100ul #C56322-AF750 100ul #C56322-Biotin 100ul #C56322-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	GRK3 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	GRK3 Antibody detects endogenous levels of total GRK3
Immunogen Description	A synthesized peptide derived from human GRK3
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	ADRBK2; BARK2; GRK3;
Accession No.	Uniprot:P35626
Calculated MW	80kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Specifically phosphorylates the agonist-occupied form of the beta-adrenergic and closely related receptors.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.