

CD52 Conjugated Antibody

Catalog No: #C56323

Package Size: #C56323-AF350 100ul #C56323-AF405 100ul #C56323-AF488 100ul #C56323-AF555 100ul #C56323-AF594 100ul #C56323-AF647 100ul #C56323-AF680 100ul #C56323-AF750 100ul #C56323-Biotin 100ul #C56323-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	CD52 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Mouse
Specificity	CD52 Antibody detects endogenous levels of total CD52
Immunogen Description	A synthesized peptide derived from human CD52
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CAMPATH 1; CD52; CDw52; He5;
Accession No.	Uniprot:Q64389
Calculated MW	8kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

May play a role in carrying and orienting carbohydrate, as well as having a more specific role.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.