

Heparanase 1 Conjugated Antibody

Catalog No: #C56324

Package Size: #C56324-AF350 100ul #C56324-AF405 100ul #C56324-AF488 100ul #C56324-AF555 100ul #C56324-AF594 100ul #C56324-AF647 100ul #C56324-AF680 100ul #C56324-AF750 100ul #C56324-Biotin 100ul #C56324-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Heparanase 1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	Heparanase 1 Antibody detects endogenous levels of total Heparanase 1
Immunogen Description	A synthesized peptide derived from human Heparanase 1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	HEP; Heparanase; Heparanase1; Hpa 1; HPA; Hpa1 HPR1; HPSE1; HSE1;
Accession No.	Uniprot:Q9Y251
Calculated MW	70kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Endoglycosidase which is a cell surface and extracellular matrix-degrading enzyme. Cleaves heparan sulfate proteoglycans (HSPGs) into heparan sulfate side chains and core proteoglycans.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.