

Alpha1 adaptin Conjugated Antibody

Catalog No: #C56328

Package Size: #C56328-AF350 100ul #C56328-AF405 100ul #C56328-AF488 100ul #C56328-AF555 100ul #C56328-AF594 100ul #C56328-AF647 100ul #C56328-AF680 100ul #C56328-AF750 100ul #C56328-Biotin 100ul #C56328-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Alpha1 adaptin Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	Alpha1-adaptin Antibody detects endogenous levels of total Alpha1-adaptin
Immunogen Description	A synthesized peptide derived from human Alpha1-adaptin
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	ADTAA; CLAPA1;
Accession No.	Uniprot:O95782
Calculated MW	108kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Component of the adaptor protein complex 2 (AP-2). Adaptor protein complexes function in protein transport via transport vesicles in different membrane traffic pathways. Adaptor protein complexes are vesicle coat components and appear to be involved in cargo selection and vesicle formation.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.