

Cytochrome P450 1A2 Conjugated Antibody

Catalog No: #C56329

Package Size: #C56329-AF350 100ul #C56329-AF405 100ul #C56329-AF488 100ul #C56329-AF555 100ul #C56329-AF594 100ul #C56329-AF647 100ul #C56329-AF680 100ul #C56329-AF750 100ul #C56329-Biotin 100ul #C56329-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Cytochrome P450 1A2 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human
Specificity	Cytochrome P450 1A2 Antibody detects endogenous levels of total Cytochrome P450 1A2
Immunogen Description	A synthesized peptide derived from human Cytochrome P450 1A2
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CP12; CYP1A2; CYP1A2; P3 450; P450 4; P450 P3;
Accession No.	Uniprot:P05177
Calculated MW	58kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Cytochromes P450 are a group of heme-thiolate monooxygenases. In liver microsomes, this enzyme is involved in an NADPH-dependent electron transport pathway.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.