

# Thyroglobulin Conjugated Antibody

Catalog No: #C56331

Package Size: #C56331-AF350 100ul #C56331-AF405 100ul #C56331-AF488 100ul #C56331-AF555 100ul #C56331-AF594 100ul #C56331-AF647 100ul #C56331-AF680 100ul #C56331-AF750 100ul #C56331-Biotin 100ul #C56331-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Thyroglobulin Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	Thyroglobulin Antibody detects endogenous levels of total Thyroglobulin
Immunogen Description	A synthesized peptide derived from human Thyroglobulin
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	AITD3; hTG; TDH3; Tg; TGN; Thyroglobulin;
Accession No.	Uniprot:P01266
Calculated MW	305kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

## Product Description

Precursor of the iodinated thyroid hormones thyroxine (T4) and triiodothyronine (T3).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.