

PEDF Conjugated Antibody

Catalog No: #C56332

Package Size: #C56332-AF350 100ul #C56332-AF405 100ul #C56332-AF488 100ul #C56332-AF555 100ul #C56332-AF594 100ul #C56332-AF647 100ul #C56332-AF680 100ul #C56332-AF750 100ul #C56332-Biotin 100ul #C56332-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	PEDF Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human
Specificity	PEDF Antibody detects endogenous levels of total PEDF
Immunogen Description	A synthesized peptide derived from human PEDF
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	EPC1; OI12; OI6; PEDF; PIG35; Serpinf1;
Accession No.	Uniprot:P36955
Calculated MW	46kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Neurotrophic protein; induces extensive neuronal differentiation in retinoblastoma cells. Potent inhibitor of angiogenesis. As it does not undergo the S (stressed) to R (relaxed) conformational transition characteristic of active serpins, it exhibits no serine protease inhibitory activity.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.