

# IGFBP7 Conjugated Antibody

Catalog No: #C56335

Package Size: #C56335-AF350 100ul #C56335-AF405 100ul #C56335-AF488 100ul #C56335-AF555 100ul #C56335-AF594 100ul #C56335-AF647 100ul #C56335-AF680 100ul #C56335-AF750 100ul #C56335-Biotin 100ul #C56335-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	IGFBP7 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	IGFBP7 Antibody detects endogenous levels of total IGFBP7
Immunogen Description	A synthesized peptide derived from human IGFBP7
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	AGM; PSF; TAF; FSTL2; IBP-7; MAC25; IGFBP-7; RAMSVPS; IGFBP-7v; IGFBPRP1;
Accession No.	Uniprot:Q16270
Calculated MW	29kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

## Product Description

Binds IGF-I and IGF-II with a relatively low affinity. Stimulates prostacyclin (PGI<sub>2</sub>) production. Stimulates cell adhesion.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.