

# Netrin 1 Conjugated Antibody

Catalog No: #C56336

Package Size: #C56336-AF350 100ul #C56336-AF405 100ul #C56336-AF488 100ul #C56336-AF555 100ul #C56336-AF594 100ul #C56336-AF647 100ul #C56336-AF680 100ul #C56336-AF750 100ul #C56336-Biotin 100ul #C56336-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Netrin 1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	Netrin 1 Antibody detects endogenous levels of total Netrin 1
Immunogen Description	A synthesized peptide derived from human Netrin 1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Netrin1; NTN1; NTN1L; Unc6;
Accession No.	Uniprot:O95631
Calculated MW	75kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

## Product Description

Netrins control guidance of CNS commissural axons and peripheral motor axons. Its association with either DCC or some UNC5 receptors will lead to axon attraction or repulsion, respectively.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.