

Loricrin Conjugated Antibody

Catalog No: #C56338

Package Size: #C56338-AF350 100ul #C56338-AF405 100ul #C56338-AF488 100ul #C56338-AF555 100ul #C56338-AF594 100ul #C56338-AF647 100ul #C56338-AF680 100ul #C56338-AF750 100ul #C56338-Biotin 100ul #C56338-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Loricrin Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	Loricrin Antibody detects endogenous levels of total Loricrin
Immunogen Description	A synthesized peptide derived from human Loricrin
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	LOR; Loricrin; LRN;
Accession No.	Uniprot:P23490
Calculated MW	38kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Major keratinocyte cell envelope protein.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.