

# Claudin 11 Conjugated Antibody

Catalog No: #C56339

Package Size: #C56339-AF350 100ul #C56339-AF405 100ul #C56339-AF488 100ul #C56339-AF555 100ul #C56339-AF594 100ul #C56339-AF647 100ul #C56339-AF680 100ul #C56339-AF750 100ul #C56339-Biotin 100ul #C56339-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Claudin 11 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	Oligodendrocyte Specific Protein Antibody detects endogenous levels of total Oligodendrocyte Specific Protein
Immunogen Description	A synthesized peptide derived from human Oligodendrocyte Specific Protein
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Claudin 11; CLD11_HUMAN; Cldn11; OSP; OTM;
Accession No.	Uniprot:O75508
Calculated MW	22kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

## Product Description

Plays a major role in tight junction-specific obliteration of the intercellular space, through calcium-independent cell-adhesion activity.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.