

# TBX21 Conjugated Antibody

Catalog No: #C56342

Package Size: #C56342-AF350 100ul #C56342-AF405 100ul #C56342-AF488 100ul #C56342-AF555 100ul #C56342-AF594 100ul #C56342-AF647 100ul #C56342-AF680 100ul #C56342-AF750 100ul #C56342-Biotin 100ul #C56342-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	TBX21 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	IF
Species Reactivity	Human
Specificity	T-bet / Tbx21 Antibody detects endogenous levels of total T-bet / Tbx21
Immunogen Description	A synthesized peptide derived from human T-bet / Tbx21
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	T bet; T box 21; T PET; TBET; TBLYM; Tbx21; TPET;
Accession No.	Uniprot:Q9UL17
Calculated MW	kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

IF:1:50-1:200

## Product Description

Transcription factor that controls the expression of the TH1 cytokine, interferon-gamma. Initiates TH1 lineage development from naive TH precursor cells both by activating TH1 genetic programs and by repressing the opposing TH2 programs.

**Note:** This product is for in vitro research use only and is not intended for use in humans or animals.