

# Mesothelin Conjugated Antibody

Catalog No: #C56344

Package Size: #C56344-AF350 100ul #C56344-AF405 100ul #C56344-AF488 100ul #C56344-AF555 100ul #C56344-AF594 100ul #C56344-AF647 100ul #C56344-AF680 100ul #C56344-AF750 100ul #C56344-Biotin 100ul #C56344-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mesothelin Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	Mesothelin Antibody detects endogenous levels of total Mesothelin
Immunogen Description	A synthesized peptide derived from human Mesothelin
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CAK1; Mesothelin; MPF; Msln; SMR; SMRP;
Accession No.	Uniprot:Q13421
Calculated MW	69,40kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

## Product Description

Membrane-anchored forms may play a role in cellular adhesion.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.