

CD48 Conjugated Antibody

Catalog No: #C56347



Package Size: #C56347-AF350 100ul #C56347-AF405 100ul #C56347-AF488 100ul #C56347-AF555 100ul #C56347-AF594 100ul #C56347-AF647 100ul #C56347-AF680 100ul #C56347-AF750 100ul #C56347-Biotin 100ul #C56347-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	CD48 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	CD48 Antibody detects endogenous levels of total CD48
Immunogen Description	A synthesized peptide derived from human CD48
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	BCM1; BLAST; BLAST1; CD48; hCD48; mCD48; MEM102; SLAMF2;
Accession No.	Uniprot:P09326
Calculated MW	43kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Ligand for CD2. Might facilitate interaction between activated lymphocytes. Probably involved in regulating T-cell activation.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.