

# Mammaglobin A Conjugated Antibody

Catalog No: #C56348

Package Size: #C56348-AF350 100ul #C56348-AF405 100ul #C56348-AF488 100ul #C56348-AF555 100ul #C56348-AF594 100ul #C56348-AF647 100ul #C56348-AF680 100ul #C56348-AF750 100ul #C56348-Biotin 100ul #C56348-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mammaglobin A Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Rat
Specificity	Mammaglobin A Antibody detects endogenous levels of total Mammaglobin A
Immunogen Description	A synthesized peptide derived from human Mammaglobin A
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Mammaglobin 1; Mammaglobin-A; MGB1; Scgb2a2; SG2A2; UGB2;
Accession No.	Uniprot:Q13296
Calculated MW	18kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

## Product Description

Belongs to the secretoglobin family. Lipophilin subfamily.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.