

IGHD Conjugated Antibody

Catalog No: #C56351

Package Size: #C56351-AF350 100ul #C56351-AF405 100ul #C56351-AF488 100ul #C56351-AF555 100ul #C56351-AF594 100ul #C56351-AF647 100ul #C56351-AF680 100ul #C56351-AF750 100ul #C56351-Biotin 100ul #C56351-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	IGHD Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	IGHD Antibody detects endogenous levels of total IGHG
Immunogen Description	A synthesized peptide derived from human IGHG
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	IGHD; IgD; Immunoglobulin delta;
Accession No.	Uniprot:P01880
Calculated MW	84kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

There are two named isoforms of IGHG produced by alternative splicing, a secreted form and a single pass type I membrane protein.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.