

TRAP Conjugated Antibody

Catalog No: #C56352

Package Size: #C56352-AF350 100ul #C56352-AF405 100ul #C56352-AF488 100ul #C56352-AF555 100ul #C56352-AF594 100ul #C56352-AF647 100ul #C56352-AF680 100ul #C56352-AF750 100ul #C56352-Biotin 100ul #C56352-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	TRAP Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	Tartrate Resistant Acid Phosphatase Antibody detects endogenous levels of total Tartrate Resistant Acid Phosphatase
Immunogen Description	A synthesized peptide derived from human Tartrate Resistant Acid Phosphatase
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	TRAP; SPENCDI; TrATPase;
Accession No.	Uniprot:P13686
Calculated MW	35,42kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Involved in osteopontin/bone sialoprotein dephosphorylation. Its expression seems to increase in certain pathological states such as Gaucher and Hodgkin diseases, the hairy cell, the B-cell, and the T-cell leukemias.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.