

Napsin A Conjugated Antibody

Catalog No: #C56353



Package Size: #C56353-AF350 100ul #C56353-AF405 100ul #C56353-AF488 100ul #C56353-AF555 100ul #C56353-AF594 100ul #C56353-AF647 100ul #C56353-AF680 100ul #C56353-AF750 100ul #C56353-Biotin 100ul #C56353-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Napsin A Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	NAPSIN A Antibody detects endogenous levels of total NAPSIN A
Immunogen Description	A synthesized peptide derived from human NAPSIN A
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	ASP4; KAP; Kdap; NAP1; NAPA; Napsa; Napsin-1; Napsin-A; Pronapsin A; SNAPA; TA01/TA02;
Accession No.	Uniprot:O96009
Calculated MW	45kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

May be involved in processing of pneumocyte surfactant precursors.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.