

DMRT1 Conjugated Antibody

Catalog No: #C56354

Package Size: #C56354-AF350 100ul #C56354-AF405 100ul #C56354-AF488 100ul #C56354-AF555 100ul #C56354-AF594 100ul #C56354-AF647 100ul #C56354-AF680 100ul #C56354-AF750 100ul #C56354-Biotin 100ul #C56354-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	DMRT1 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	DMRT1 Antibody detects endogenous levels of total DMRT1
Immunogen Description	A synthesized peptide derived from human DMRT1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CT154; dmt1; DMT1;
Accession No.	Uniprot:Q9Y5R6
Calculated MW	34kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

May be required for testis development.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.