

MSMB Conjugated Antibody

Catalog No: #C56355



Package Size: #C56355-AF350 100ul #C56355-AF405 100ul #C56355-AF488 100ul #C56355-AF555 100ul #C56355-AF594 100ul #C56355-AF647 100ul #C56355-AF680 100ul #C56355-AF750 100ul #C56355-Biotin 100ul #C56355-Conjugated 50ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MSMB Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	Prostate Secretory Protein/PSP Antibody detects endogenous levels of total Prostate Secretory Protein/PSP
Immunogen Description	A synthesized peptide derived from human Prostate Secretory Protein/PSP
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	MSP; PSP; IGBF; MSPB; PN44; PRPS; HPC13; PSP57; PSP94; PSP-94;
Accession No.	Uniprot:P08118
Calculated MW	15kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Specific receptors for this protein are found on spermatozoa and in the prostate.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.