

K Cadherin Conjugated Antibody

Catalog No: #C56356



Package Size: #C56356-AF350 100ul #C56356-AF405 100ul #C56356-AF488 100ul #C56356-AF555 100ul #C56356-AF594 100ul #C56356-AF647 100ul #C56356-AF680 100ul #C56356-AF750 100ul #C56356-Biotin 100ul #C56356-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	K Cadherin Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	K Cadherin Antibody detects endogenous levels of total K Cadherin
Immunogen Description	A synthesized peptide derived from human K Cadherin
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CAD6; Cdh6; K-cadherin; KCAD; Kidney cadherin;
Accession No.	Uniprot:P55285
Calculated MW	130kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Cadherins are calcium dependent cell adhesion proteins. They preferentially interact with themselves in a homophilic manner in connecting cells; cadherins may thus contribute to the sorting of heterogeneous cell types.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.