

TSH beta Conjugated Antibody

Catalog No: #C56357

Package Size: #C56357-AF350 100ul #C56357-AF405 100ul #C56357-AF488 100ul #C56357-AF555 100ul #C56357-AF594 100ul #C56357-AF647 100ul #C56357-AF680 100ul #C56357-AF750 100ul #C56357-Biotin 100ul #C56357-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	TSH beta Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	TSH beta Antibody detects endogenous levels of total TSH beta
Immunogen Description	A synthesized peptide derived from human TSH beta
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CHNG4; TSHB;
Accession No.	Uniprot:P01222
Calculated MW	16kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Various hormones are secreted from the anterior pituitary during development and growth. Thyroid Stimulating Hormone (TSH) is secreted from the anterior pituitary and controls thyroid structure and metabolism

Note: This product is for in vitro research use only and is not intended for use in humans or animals.