

# IMP3 Conjugated Antibody

Catalog No: #C56358

Package Size: #C56358-AF350 100ul #C56358-AF405 100ul #C56358-AF488 100ul #C56358-AF555 100ul #C56358-AF594 100ul #C56358-AF647 100ul #C56358-AF680 100ul #C56358-AF750 100ul #C56358-Biotin 100ul #C56358-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	IMP3 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	IMP3 Antibody detects endogenous levels of total IMP3
Immunogen Description	A synthesized peptide derived from human IMP3
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CT98; hKOC; IGF2BP3; IMP 3; KOC1; VICKZ 3;
Accession No.	Uniprot:O00425
Calculated MW	69kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

## Product Description

RNA-binding protein that act as a regulator of mRNA translation and stability. Binds to the 5'-UTR of the insulin-like growth factor 2 (IGF2) mRNAs. Binds to sequences in the 3'-UTR of CD44 mRNA.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.