

PGDH Conjugated Antibody

Catalog No: #C56360

Package Size: #C56360-AF350 100ul #C56360-AF405 100ul #C56360-AF488 100ul #C56360-AF555 100ul #C56360-AF594 100ul #C56360-AF647 100ul #C56360-AF680 100ul #C56360-AF750 100ul #C56360-Biotin 100ul #C56360-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	PGDH Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human
Specificity	Prostaglandin dehydrogenase 1 Antibody detects endogenous levels of total Prostaglandin dehydrogenase 1
Immunogen Description	A synthesized peptide derived from human Prostaglandin dehydrogenase 1
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	15-PGDH; Hpgd; PGDH; PGDH1; PHOAR1; SDR36C1;
Accession No.	Uniprot:P15428
Calculated MW	29kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Prostaglandin inactivation. Contributes to the regulation of events that are under the control of prostaglandin levels. Catalyzes the NAD-dependent dehydrogenation of lipoxin A4 to form 15-oxo-lipoxin A4. Inhibits in vivo proliferation of colon cancer cells.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.