

Uroplakin3a Conjugated Antibody

Catalog No: #C56361

Package Size: #C56361-AF350 100ul #C56361-AF405 100ul #C56361-AF488 100ul #C56361-AF555 100ul #C56361-AF594 100ul #C56361-AF647 100ul #C56361-AF680 100ul #C56361-AF750 100ul #C56361-Biotin 100ul #C56361-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Uroplakin3a Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	IF, FC
Species Reactivity	Human
Specificity	Uroplakin III Antibody detects endogenous levels of total Uroplakin III
Immunogen Description	A synthesized peptide derived from human Uroplakin III
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	UP3a; UPIII; UPIIIA; UPK3; UPK3A; Uroplakin3A; UroplakinIII;
Accession No.	Uniprot:O75631
Calculated MW	kDa
Storage	Store at 4°C in dark for 6 months

Application Details

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Component of the asymmetric unit membrane (AUM); a highly specialized biomembrane elaborated by terminally differentiated urothelial cells. May play an important role in AUM-cytoskeleton interaction in terminally differentiated urothelial cells.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.