

Cathepsin S Conjugated Antibody

Catalog No: #C56363

Package Size: #C56363-AF350 100ul #C56363-AF405 100ul #C56363-AF488 100ul #C56363-AF555 100ul #C56363-AF594 100ul #C56363-AF647 100ul #C56363-AF680 100ul #C56363-AF750 100ul #C56363-Biotin 100ul #C56363-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Cathepsin S Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	Cathepsin S Antibody detects endogenous levels of total Cathepsin S
Immunogen Description	A synthesized peptide derived from human Cathepsin S
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CTSS; Cathepsin S; Cat-s; CATS;
Accession No.	Uniprot:P25774
Calculated MW	25kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Thiol protease. Key protease responsible for the removal of the invariant chain from MHC class II molecules. The bond-specificity of this proteinase is in part similar to the specificities of cathepsin L and cathepsin N.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.