

PMP70 Conjugated Antibody

Catalog No: #C56364

Package Size: #C56364-AF350 100ul #C56364-AF405 100ul #C56364-AF488 100ul #C56364-AF555 100ul #C56364-AF594 100ul #C56364-AF647 100ul #C56364-AF680 100ul #C56364-AF750 100ul #C56364-Biotin 100ul #C56364-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	PMP70 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	PMP70 Antibody detects endogenous levels of total PMP70
Immunogen Description	A synthesized peptide derived from human PMP70
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	ABCD3; ABC43; PMP70; PXMP1; ZWS2;
Accession No.	Uniprot:P28288
Calculated MW	70kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Probable transporter. The nucleotide-binding fold acts as an ATP-binding subunit with ATPase activity.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.