

NPRC Conjugated Antibody

Catalog No: #C56370

Package Size: #C56370-AF350 100ul #C56370-AF405 100ul #C56370-AF488 100ul #C56370-AF555 100ul #C56370-AF594 100ul #C56370-AF647 100ul #C56370-AF680 100ul #C56370-AF750 100ul #C56370-Biotin 100ul #C56370-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	NPRC Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF, FC
Species Reactivity	Human Mouse Rat
Specificity	Natriuretic Peptide Receptor C Antibody detects endogenous levels of total Natriuretic Peptide Receptor C
Immunogen Description	A synthesized peptide derived from human Natriuretic Peptide Receptor C
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	ANPC; ANPRC; GUCY2B; NPR 3; NPR C; Npr3; NPRC;
Accession No.	Uniprot:P17342
Calculated MW	70kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

FC: 1:50-1:200

Product Description

Receptor for the natriuretic peptide hormones, binding with similar affinities atrial natriuretic peptide NPPA/ANP, brain natriuretic peptide NPPB/BNP, and C-type natriuretic peptide NPPC/CNP.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.