

## CRALBP Conjugated Antibody

Catalog No: #C56371



Package Size: #C56371-AF350 100ul #C56371-AF405 100ul #C56371-AF488 100ul #C56371-AF555 100ul #C56371-AF594 100ul #C56371-AF647 100ul #C56371-AF680 100ul #C56371-AF750 100ul #C56371-Biotin 100ul #C56371-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CRALBP Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Mouse Rat
Specificity	CRALBP Antibody detects endogenous levels of total CRALBP
Immunogen Description	A synthesized peptide derived from human CRALBP
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	RLBP1;
Accession No.	Uniprot:P12271
Calculated MW	36kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

## Product Description

Participates in the regeneration of active 11-cis-retinol and 11-cis-retinaldehyde, from the inactive 11-trans products of the rhodopsin photocycle and in the de novo synthesis of these retinoids from 11-trans metabolic precursors.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.