

# OLFM4 Conjugated Antibody

Catalog No: #C56375

Package Size: #C56375-AF350 100ul #C56375-AF405 100ul #C56375-AF488 100ul #C56375-AF555 100ul #C56375-AF594 100ul #C56375-AF647 100ul #C56375-AF680 100ul #C56375-AF750 100ul #C56375-Biotin 100ul #C56375-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	OLFM4 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human Mouse Rat
Specificity	OLFM4 Antibody detects endogenous levels of total OLFM4
Immunogen Description	A synthesized peptide derived from human OLFM4
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	GC1; GW112; hGC 1; hOLFID; olfactomedin 4; OlfD; OLFM4; OLM4;
Accession No.	Uniprot:Q6UX06
Calculated MW	54kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF:1:50-1:200

## Product Description

An antiapoptotic factor that promotes tumor growth. Promotes proliferation of pancreatic cancer cells by favoring the transition from the S to G2/M phase. Facilitates cell adhesion.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.