

Alkaline phosphatase Conjugated Antibody

Catalog No: #C56377



Package Size: #C56377-AF350 100ul #C56377-AF405 100ul #C56377-AF488 100ul #C56377-AF555 100ul #C56377-AF594 100ul
 #C56377-AF647 100ul #C56377-AF680 100ul #C56377-AF750 100ul #C56377-Biotin 100ul #C56377-Conjugated 50ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Alkaline phosphatase Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	intestinal alkaline phosphatase Antibody detects endogenous levels of total intestinal alkaline phosphatase
Immunogen Description	A synthesized peptide derived from human intestinal alkaline phosphatase
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Alpi; Glycerophosphatase; IAP;
Accession No.	Uniprot:P09923
Calculated MW	68kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

In most mammals there are four different isozymes: placental, placental-like, intestinal and tissue non-specific (liver/bone/kidney).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.