

# Cytokeratin 6 Conjugated Antibody

Catalog No: #C56381

Package Size: #C56381-AF350 100ul #C56381-AF405 100ul #C56381-AF488 100ul #C56381-AF555 100ul #C56381-AF594 100ul #C56381-AF647 100ul #C56381-AF680 100ul #C56381-AF750 100ul #C56381-Biotin 100ul #C56381-Conjugated 50ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Cytokeratin 6 Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	CK6 Antibody detects endogenous levels of total CK6
Immunogen Description	A synthesized peptide derived from human CK6
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CK 6A; CK 6B; CK 6C; CK 6D; CK 6E; CK-6C; CK-6E; K6a keratin; K6b keratin; K6C; K6c keratin; K6d keratin; K6e keratin; Keratin; KRT6A; KRT6B; KRT6C; KRT6D; KRT6E;
Accession No.	Uniprot:P02538/P04259/P48668
Calculated MW	60kDa
Storage	Store at 4°C in dark for 6 months

## Application Details

WB: 1:50-1:200

IF: 1:50-1:200

## Product Description

There are at least six isoforms of human type II keratin-6 (K6). There are two types of cytoskeletal and microfibrillar keratin, I (acidic) and II (neutral to basic) (40-55 and 56-70 kDa, respectively).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.