

TPO Conjugated Antibody

Catalog No: #C56389

Package Size: #C56389-AF350 100ul #C56389-AF405 100ul #C56389-AF488 100ul #C56389-AF555 100ul #C56389-AF594 100ul #C56389-AF647 100ul #C56389-AF680 100ul #C56389-AF750 100ul #C56389-Biotin 100ul #C56389-Conjugated 50ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	TPO Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB, IF
Species Reactivity	Human
Specificity	TPO Antibody detects endogenous levels of total TPO
Immunogen Description	A synthesized peptide derived from human TPO
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	MSA; TDH2A; TPO; TPX;
Accession No.	Uniprot:P07202
Calculated MW	103kDa
Storage	Store at 4°C in dark for 6 months

Application Details

WB: 1:50-1:200

IF:1:50-1:200

Product Description

Iodination and coupling of the homonogenic tyrosines in thyroglobulin to yield the thyroid hormones T(3) and T(4).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.